
Automation Solutions Products

General Safety Guide 

Contents

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1 English

This guide describes the general safety precautions, intended product use, and the list of safety labels for the Automation Solutions products. The user documentation supplied with your product might provide additional specific safety requirements for the product.

This guide contains the following topics:

- “Before using the Automation Solutions product” on page 2
- “Third-party products” on page 2
- “Safety notices” on page 2
- “Intended product use” on page 2
- “Safety labels” on page 3
- “Ergonomics” on page 4
- “Electrical hazards” on page 4
- “Radiation hazards” on page 4
- “Chemical hazards” on page 5
- “Gas hazards” on page 6
- “Heat hazards” on page 7
- “Moving-parts hazards” on page 7

For other user information, you can search the product knowledge base as follows:

- Use the help menu in the Automation Solutions lab automation system software.
- Go to the Automation Solutions website at www.agilent.com/chem/askb.

Original Instructions

Before using the Automation Solutions product

You must be aware of the potential hazards and understand how to avoid being exposed to them. In addition, make sure you are properly trained in:

- General laboratory safety
- The correct and safe operation of the product
- The correct and safe operation of other lab automation systems or components used in combination with the product

Third-party products

Third-party products that are integrated in Automation Solutions lab automation systems might contain electrical, chemical, or other hazards. For information about safety for third-party products, see the relevant product documentation.

Safety notices

Pay special attention to the following safety notices in all product documentation.

WARNING A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

CAUTION A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

Intended product use

WARNING Do not remove the product exterior covers or otherwise disassemble the product. Doing so can expose you to hazards that could cause serious injury and damage the product.

WARNING Using controls, making adjustments, or performing procedures other than those specified in the product user documentation can expose you to hazardous voltage.

Agilent Technologies products must only be used in the manner described in the Agilent Technologies product user documentation. Any other use might result in damage to the product or personal injury. Agilent Technologies is not responsible for any damages caused, in whole or in part, by improper use of the products, unauthorized alterations, adjustments or modifications to the products, failure to comply with procedures in Agilent Technologies product user documentation, or

use of the products in violation of applicable laws, rules or regulations. Except as otherwise expressly provided in Agilent Technologies product user documentation, any alteration, adjustment, or modification to the products will void the product warranty.

The product is not intended or approved for diagnosis of disease in humans or animals. You assume full responsibility for obtaining any regulatory approvals required for such use and assume all liability in connection therewith.

Safety labels

The following table lists the common safety labels you might find on the Automation Solutions products.

Symbol	Description
	Indicates that you must read the accompanying instructions (for example, the safety guide) for more information before proceeding.
	Indicates hazardous voltages.
	Indicates pinch, crush, or cut hazard.
	
	
	Indicates laser hazard.
	Indicates hot surface hazard.
	Indicates that to prevent serious injury, guards, interlocks, and safety devices should be restored to correct operation when installing parts or making repairs. Refer to user or service documentation for additional information.
	Indicates protective conductor terminal, which is bonded to conductive parts of an equipment for protection against electric shock in case of a fault and is intended to be connected to an external protective earthing system.
	Indicates frame or chassis terminal, which is bonded to conductive parts of an equipment for safety purposes.

Symbol	Description
	Indicates that you must not discard this electrical or electronic product in domestic household waste.

Ergonomics

The product might be large and weigh more than 18 kg (40 lb).

WARNING To avoid injury to personnel and potential equipment damage, use sufficient assistance and proper lifting techniques when moving or lifting the product.

See the computer manufacturer documentation for instructions on the correct positioning of the computer, monitor, keyboard, and mouse.

Electrical hazards

Hazardous-voltage electronics can be found within the product. Under normal operating conditions, you are protected from exposure to the hazardous voltage.

WARNING Exposure to the product internal electronics can cause severe injury. Do not remove the product's exterior cover panels or try to gain access to the internal electronics.

WARNING Ensure that the power cords are in good condition and are not frayed. Use of frayed or damaged power cords can cause injury.

CAUTION Verify that the voltage range and frequency of your power distribution match the power specifications of the individual product. Connecting the product to incompatible circuits can cause damage to the product.

CAUTION Always use the supplied power cords. Use of incorrect power cords can cause damage to the product.

Hazardous-voltage electronics can also be found in the controlling computer. See the computer manufacturer documentation for the hazard warnings. Make sure you follow the instructions on the safe operation of the computer.

Radiation hazards

The Automation Solutions product can contain one or more barcode readers that contain lasers for reading barcodes. For example, the Mirrored Barcode Reader, Auxiliary Barcode Reader, and the optional barcode reader on the Microplate Labeler contain class II lasers.

WARNING Do not look directly at the laser beam. Looking directly at the laser beam can result in serious eye injury.

WARNING Do not disassemble barcode reader sensor heads. Laser emission from the reader is not automatically stopped if the sensor head is disassembled, and can cause potential injury.

WARNING Unless you are following a procedure in a Agilent Technologies user documentation, do not touch the white test button on the side of the sensor head of an auxiliary barcode reader. The white test button turns on the laser and can result in unsafe exposure to the beam.

Barcode readers emit light for up to 0.5 seconds only when taking a reading. The potential hazard is present only during protocol runs. When a microplate is being read, the beam should not pass beyond the microplate. Operating the barcode readers without microplates increases the risk of exposure.

Third-party products that are integrated in Automation Solutions lab automation systems might also contain lasers, such as Zeiss Multimode readers and Liconic incubators. For information about laser safety for third-party products, see the relevant product documentation.

Chemical hazards

General guidelines

Some chemicals used when working with the Automation Solutions product can be hazardous. Make sure you:

- Follow standard laboratory procedures and cautions when working with chemicals.
- Follow safety directions in the Material Safety Data Sheet (MSDS) for each chemical.
- Follow your local, state, and federal safety regulations when using and disposing of the chemicals.

Compatible chemicals

The Automation Solutions product surfaces are designed to be compatible with small volumes of aqueous solutions, common biological buffers, solvents, and common reducing agents.

The following chemicals were tested for brief exposures:

- Water (deionized, ultrapure type 1)
- Dimethyl sulfoxide (100%)
- Strong acid, such as hydrochloric acid (100 mM HCl) pH 2-3
- Strong base (HPO₄ or sodium hydroxide, 100 mM) pH 12-13
- Ethanol (95% alcohol)
- Isopropanol (isopropyl alcohol)
- Methanol (95% alcohol)
- Acetonitrile (100%)
- Chloroform (100%)
- Hydrogen peroxide (35%)
- Tetrahydrofuran
- Dichloromethane
- Toluene
- Acetone

Agilent Technologies recommends that you employ an Environmental Health and Safety professional on site to confirm the compatibility of other chemicals.

Gas hazards

High-pressure gas

Compressed air is used to move components of some Automation Solutions products.

WARNING Working with open, charged air lines can result in injury. Turn off the compressed air line when disconnecting or reconnecting devices that use compressed air. Contact your facilities department or Automation Solutions Technical Support with questions about setting up the air line.

Follow the local, state, and federal safety codes for the placement and mounting of gas cylinders. For example, you might have to attach a standard cylinder bracket to a solid permanent structure to meet or exceed all local seismic and safety requirements.

Always use good laboratory practices when handling high-pressure cylinders. Make sure you follow any instructions provided with the cylinders.

Low-oxygen hazards

Some Automation Solutions products use inert gas to displace air. Under normal operating conditions you are protected from hazardous exposure to the gas. However, you must be aware of the hazards and how to avoid them.

WARNING Inert gases such as nitrogen and argon are odorless, colorless, and nontoxic gases that can cause rapid suffocation by displacing air. Use the inert gas only in a well-ventilated lab area. Turn off the inert gas sources when the product is not in use.

Heat hazards

Some Automation Solutions products become hot during operation.

WARNING Do not touch the product while it is in operation. The hot surfaces can cause burn injury. Allow the device to cool before servicing.

Moving-parts hazards

Many Automation Solutions products have components that move during operation. See the product user documentation for details on the axes of motion and the reach of the moving parts.


WARNING The product's moving parts can injure you if you deviate from the procedures in the product user documentation. To avoid potential injury, keep clear of the product while it is in motion.

Many Automation Solutions products have safety features, such as physical shields, the Light Curtain, and safety interlock systems, that are designed to protect operators from moving-parts hazards. It is the responsibility of every operator to use the safety features and follow warnings and safety labels. If you have questions about the safety features, contact Automation Solutions Technical Support.

WARNING Changing or modifying the product safety features might prevent the safe operation of the product, invalidate its safety compliance, and lead to personal injury or property damage. Any customer who does not use the supplied safety equipment or who modifies the supplied safety equipment assumes full responsibility for providing an appropriate level of safety for its operators and for providing the applicable safety compliance marking and documentation.

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© Agilent Technologies, Inc. 2014
Agilent Technologies Inc.